## Battery Care, Maintenance, and Safety Guidelines

### Introduction

Maintaining a healthy battery is essential for both your safety and the longevity of your devices. This document provides guidelines on how to care for and store batteries, recognize when they are nearing end-of-life, and properly handle and dispose of swollen or otherwise compromised batteries. Following these best practices can extend battery life, optimize your device's performance, and protect yourself and others from potential hazards.

### Battery basics

- Charge to 80-90% regularly; 100% isn't necessary.
- Avoid running the battery fully flat (0%); recharge at a minimum of 40%.
- Keep your device and battery out of extreme heat or cold.

### Product Storage

- If a product is not used, it should be kept at 40% to 80%. This will prevent battery degradation and increase storage safety.
- Store product in a controlled environment between 15-25°C / 60-80°F
- The battery level should be checked regularly, and the product should be recharged if the battery level is under 40%.

### End-of-life (EOL) signs for batteries

- The device shuts down due to low battery, even if it should have plenty of time left.
- Battery charge drops randomly or rapidly.
- Charging finishes quickly when you know it should have taken longer.
- Battery capacity drops rapidly with little or no warning.
- Discharge rates are very rapid or uneven.
- The battery pack gets very hot during charging or discharging.
- The battery begins to show deformity, bulging, or swelling.

### How to handle a battery that is starting to swell or that is already swollen

1. **Do not** press on the swollen battery or its casing.
2. **Stop charging** a device that has a swollen battery.
3. Power down the device and **do not** attempt to use it.
4. **Do not** remove the battery if the device feels hot or smells.
5. **Do not** use tools of any type to pry on or against the battery when removing.
6. If your device has a removable battery, take pictures of it and the device, if possible, and discard them safely.
7. Order a new battery from HumanWare. **Do not** purchase Humanware devices or batteries from unauthorized resellers. Use of these will void any warranty, written or implied.
8. Send pictures of the battery or device to HumanWare customer support. Include the product model, serial number, approximate purchase date, and battery replacement date if it was ever replaced.

Here are more basics on [what to do with a swollen battery from iFixit](#).

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